

Recombinant Human EPHB6

Catalog No: #GP11166



Package Size: #GP11166-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EPHB6
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human EPH receptor B6
Target Name	EPH receptor B6
Other Names	HEP
Accession No.	Swissprot:O15197Gene Accession:BC110606
Uniprot	O15197
GeneID	2051;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of a family of transmembrane proteins that function as receptors for ephrin-B family proteins. Unlike other members of this family, the encoded protein does not contain a functional kinase domain. Activity of this protein can influence cell adhesion and migration. Expression of this gene is downregulated during tumor progression, suggesting that the protein may suppress tumor invasion and metastasis. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only