

Recombinant Human NEK11

Catalog No: #GP11181



Package Size: #GP11181-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NEK11
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-466 amino acids of human NIMA-related kinase 11
Target Name	NIMA-related kinase 11
Accession No.	Swissprot:Q8NG66Gene Accession:BC028587
Uniprot	Q8NG66
GeneID	79858;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the never in mitosis gene A family of kinases. The encoded protein localizes to the nucleoli, and may function with NEK2A in the S-phase checkpoint. The encoded protein appears to play roles in DNA replication and response to genotoxic stress. Alternatively spliced transcript variants have been described.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only