

Recombinant Human PACSIN2

Catalog No: #GP11184



Package Size: #GP11184-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PACSIN2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human protein kinase C and casein kinase substrate in neurons 2
Target Name	protein kinase C and casein kinase substrate in neurons 2
Other Names	SDPII
Accession No.	Swissprot:Q9UNF0Gene Accession:BC008037
Uniprot	Q9UNF0
GeneID	11252;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of the protein kinase C and casein kinase substrate in neurons family. The encoded protein is involved in linking the actin cytoskeleton with vesicle formation by regulating tubulin polymerization. Alternative splicing results in multiple transcript variants.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only