

Recombinant Human HK3

Catalog No: #GP11188



Package Size: #GP11188-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HK3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 268-515 amino acids of human hexokinase 3 (white cell)
Target Name	hexokinase 3 (white cell)
Other Names	HXK3; HKIII
Accession No.	Swissprot:P52790Gene Accession:BC028129
Uniprot	P52790
GeneID	3101;
Storage	-20~-80°C, pH 7.6 PBS

Background

Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. This gene encodes hexokinase 3. Similar to hexokinases 1 and 2, this allosteric enzyme is inhibited by its product glucose-6-phosphate.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only