

Recombinant Human CBY1

Catalog No: #GP11199



Package Size: #GP11199-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CBY1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	chibby homolog 1 (Drosophila)
Other Names	CBY; arb1; PGEA1; C22orf2; PIGEA14; PIGEA-14; HS508I15A
Accession No.	Swissprot:Q9Y3M2Gene Accession:BC016139
Uniprot	Q9Y3M2
GeneID	25776;
Storage	-20~-80°C, pH 7.6 PBS

Background

Beta-catenin is a transcriptional activator and oncoprotein involved in the development of several cancers. The protein encoded by this gene interacts directly with the C-terminal region of beta-catenin, inhibiting oncogenic beta-catenin-mediated transcriptional activation by competing with transcription factors for binding to beta-catenin. Two transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only