

Recombinant Human NT5E

Catalog No: #GP11208



Package Size: #GP11208-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NT5E
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-252 amino acids of human 5'-nucleotidase, ecto (CD73)
Target Name	5'-nucleotidase, ecto (CD73)
Other Names	NT; eN; NT5; NTE; eNT; CD73; E5NT; CALJA
Accession No.	Swissprot:P21589Gene Accession:BC015940
Uniprot	P21589
GenID	4907;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a plasma membrane protein that catalyzes the conversion of extracellular nucleotides to membrane-permeable nucleosides. The encoded protein is used as a determinant of lymphocyte differentiation. Defects in this gene can lead to the calcification of joints and arteries. Two transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only