

Recombinant Human XPO1

Catalog No: #GP11212



Package Size: #GP11212-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human XPO1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human exportin 1 (CRM1 homolog, yeast)
Target Name	exportin 1 (CRM1 homolog, yeast)
Other Names	emb; CRM1; exp1
Accession No.	Swissprot:O14980Gene Accession:BC032847
Uniprot	O14980
GeneID	7514;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene mediates leucine-rich nuclear export signal (NES)-dependent protein transport. Exportin 1 specifically inhibits the nuclear export of Rev and U snRNAs. It is involved in the control of several cellular processes by controlling the localization of cyclin B, MPAK, and MAPKAP kinase 2. This protein also regulates NFAT and AP-1.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only