

Recombinant Human MTMR4

Catalog No: #GP11213



Package Size: #GP11213-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MTMR4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 153-400 amino acids of human myotubularin related protein 4
Target Name	myotubularin related protein 4
Other Names	ZFYVE11; FYVE-DSP2
Accession No.	Swissprot:Q9NYA4Gene Accession:BC035609
Uniprot	Q9NYA4
GenID	9110;
Storage	-20~-80°C, pH 7.6 PBS

Background

MTMR4(myotubularin-related protein-4) cDNA and cDNAs encoding 3 other novel members of the myotubularin protein family. Dephosphorylates proteins phosphorylated on Ser, Thr, and Tyr residues and low molecular weight phosphatase substrate para-nitrophenylphosphate. Phosphorylates phosphatidylinositol 3,4,5-trisphosphate (PIP3).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only