

Recombinant Human KHK

Catalog No: #GP11214



Package Size: #GP11214-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KHK
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ketohehexokinase (fructokinase)
Accession No.	Swissprot:P50053Gene Accession:BC006233
Uniprot	P50053
GeneID	3795;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes ketohehexokinase that catalyzes conversion of fructose to fructose-1-phosphate. The product of this gene is the first enzyme with a specialized pathway that catabolizes dietary fructose. Alternatively spliced transcript variants encoding different isoforms have been identified.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only