

Recombinant Human PGK1

Catalog No: #GP11217



Package Size: #GP11217-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PGK1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human phosphoglycerate kinase 1
Target Name	phosphoglycerate kinase 1
Other Names	PGKA; MIG10
Accession No.	Swissprot:P00558Gene Accession:BC113568
Uniprot	P00558
GeneID	5230;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a glycolytic enzyme that catalyzes the conversion of 1,3-diphosphoglycerate to 3-phosphoglycerate. The encoded protein may also act as a cofactor for polymerase alpha. This gene lies on the X-chromosome, while a related pseudogene also has been found on the X-chromosome and another on chromosome 19.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only