

Recombinant Human ASB4

Catalog No: #GP11232



Package Size: #GP11232-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ASB4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 280 amino acids of human ankyrin repeat and SOCS box containing 4
Target Name	ankyrin repeat and SOCS box containing 4
Other Names	ASB-4
Accession No.	Swissprot:Q9Y574Gene Accession:BC113479
Uniprot	Q9Y574
GenID	51666;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. They contain ankyrin repeat sequence and SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Multiple alternatively spliced transcript variants have been described for this gene but some of the full length sequences are not known.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only