

Recombinant Human ADCY3

Catalog No: #GP11244



Package Size: #GP11244-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ADCY3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 290 amino acids of human adenylate cyclase 3
Target Name	adenylate cyclase 3
Other Names	AC3
Accession No.	Swissprot:O60266Gene Accession:BC126235
Uniprot	O60266
GeneID	109;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes adenylyl cyclase 3 which is a membrane-associated enzyme and catalyzes the formation of the secondary messenger cyclic adenosine monophosphate (cAMP). This protein appears to be widely expressed in various human tissues and may be involved in a number of physiological and pathophysiological metabolic processes.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only