

## Recombinant Human ORM2

Catalog No: #GP11245



Package Size: #GP11245-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ORM2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 19-201 amino acids of human orosomuroid 2
Target Name	orosomuroid 2
Other Names	AGP2; AGP-B; AGP-B'
Accession No.	Swissprot:P19652Gene Accession:BC015964
Uniprot	P19652
GeneID	5005;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only