

Recombinant Human STK17B

Catalog No: #GP11251



Package Size: #GP11251-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human STK17B |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 20-100 amino acids of human STK17B |
| Target Name | serine/threonine kinase 17b |
| Other Names | DRAK2 |
| Accession No. | Swissprot:O94768Gene Accession:BC016040 |
| Uniprot | O94768 |
| GeneID | 9262; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

Phosphorylates myosin light chains . Acts as a positive regulator of apoptosis.Highly expressed in placenta, lung, pancreas. Lower levels in heart, brain, liver, skeletal muscle and kidney.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only