

Recombinant Human PDS5B

Catalog No: #GP11255



Package Size: #GP11255-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PDS5B
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 280-529 amino acids of human PDS5B
Target Name	PDS5 cohesin associated factor B
Other Names	AS3; APRIN; CG008
Accession No.	Swissprot:Q9NTI5Gene Accession:BC039256
Uniprot	Q9NTI5
GeneID	23047;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that interacts with the conserved protein complex termed cohesin. The cohesin complex holds together sister chromatids and facilitates accurate chromosome segregation during mitosis and meiosis. This protein is also a negative regulator of cell proliferation and may be a tumor-suppressor gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only