

Recombinant Human ARL11

Catalog No: #GP11256



Package Size: #GP11256-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ARL11
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ADP-ribosylation factor-like 11
Other Names	ARLTS1
Accession No.	Swissprot:Q969Q4Gene Accession:BC013150
Uniprot	Q969Q4
GeneID	115761;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a tumor suppressor related to the ADP-ribosylation factor (ARF) family of proteins. The encoded protein may play a role in apoptosis in a caspase-dependent manner. Polymorphisms in this gene have been associated with some familial cancers.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only