

# Recombinant Human ASB3

Catalog No: #GP11259



Package Size: #GP11259-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ASB3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 310-455 amino acids of human ankyrin repeat and SOCS box containing 3
Target Name	ankyrin repeat and SOCS box containing 3
Other Names	ASB-3
Accession No.	Swissprot:Q9Y575Gene Accession:BC006488
Uniprot	Q9Y575
GenelD	100302652;51130;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. They contain ankyrin repeat sequence and SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Alternatively spliced transcript variants have been described for this gene.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only