

Recombinant Human ASB10

Catalog No: #GP11260



Package Size: #GP11260-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ASB10
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 310-445 amino acids of human ASB10
Target Name	ankyrin repeat and SOCS box containing 10
Other Names	GLC1F
Accession No.	Swissprot:Q8WXI3Gene Accession:BC126351
Uniprot	Q8WXI3
GeneID	136371;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. The SOCS box serves to couple suppressor of cytokine signaling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Multiple alternatively spliced transcript variants have been described for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only