

Recombinant Human TP53BP2

Catalog No: #GP11262



Package Size: #GP11262-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TP53BP2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human tumor protein p53 binding protein, 2
Target Name	tumor protein p53 binding protein, 2
Other Names	BBP; 53BP2; ASPP2; P53BP2; PPP1R13A
Accession No.	Swissprot:Q13625Gene Accession:BC040247
Uniprot	Q13625
GeneID	7159;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the ASPP (apoptosis-stimulating protein of p53) family of p53 interacting proteins. The protein contains four ankyrin repeats and an SH3 domain involved in protein-protein interactions. It is localized to the perinuclear region of the cytoplasm, and regulates apoptosis and cell growth through interactions with other regulatory molecules including members of the p53 family. Multiple transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only