

Recombinant Human SH3GLB1

Catalog No: #GP11281



Package Size: #GP11281-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SH3GLB1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 27-261 amino acids of human SH3GLB1
Target Name	SH3-domain GRB2-like endophilin B1
Other Names	Bif-1; CGI-61; PPP1R70; dJ612B15.2
Accession No.	Swissprot:Q9Y371 Gene Accession:BC007455
Uniprot	Q9Y371
GeneID	51100;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a SRC homology 3 domain-containing protein. The encoded protein interacts with the proapoptotic member of the Bcl-2 family, Bcl-2-associated X protein (Bax) and may be involved in regulating apoptotic signaling pathways. This protein may also be involved in maintaining mitochondrial morphology. Alternate splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only