

Recombinant Human BIN3

Catalog No: #GP11283



Package Size: #GP11283-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BIN3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 9-232 amino acids of human bridging integrator 3
Target Name	bridging integrator 3
Accession No.	Swissprot:Q9NQY0Gene Accession:BC009824
Uniprot	Q9NQY0
GeneID	55909;
Storage	-20~-80°C, pH 7.6 PBS

Background

The product of this gene is a member of the BAR domain protein family. The encoded protein is comprised solely of a BAR domain which is predicted to form coiled-coil structures and proposed to mediate dimerization, sense and induce membrane curvature, and bind small GTPases. BAR domain proteins have been implicated in endocytosis, intracellular transport, and a diverse set of other processes

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only