

Recombinant Human KEL

Catalog No: #GP11287



Package Size: #GP11287-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human KEL
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human Kell blood group, metallo-endopeptidase
Target Name	Kell blood group, metallo-endopeptidase
Other Names	ECE3; CD238
Accession No.	Swissprot:P23276Gene Accession:BC003135
Uniprot	P23276
GenID	3792;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a type II transmembrane glycoprotein that is the highly polymorphic Kell blood group antigen. The Kell glycoprotein links via a single disulfide bond to the XK membrane protein that carries the Kx antigen. The encoded protein contains sequence and structural similarity to members of the neprilysin (M13) family of zinc endopeptidases.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only