

Recombinant Human C3AR1

Catalog No: #GP11299



Package Size: #GP11299-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Recombinant Human C3AR1 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 161-340 amino acids of human C3AR1 |
| Target Name | complement C3a receptor 1 |
| Other Names | AZ3B; C3AR; HNFAG09 |
| Accession No. | Swissprot:Q16581Gene Accession:BC020742 |
| Uniprot | Q16581 |
| GeneID | 719; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

C3a is an anaphylatoxin released during activation of the complement system. The protein encoded by this gene is an orphan G protein-coupled receptor for C3a. Binding of C3a by the encoded receptor activates chemotaxis, granule enzyme release, superoxide anion production, and bacterial opsonization.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only