

## PRELID1 Antibody

Catalog No: #42978

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PRELID1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total PRELID1 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human PRELID1
Target Name	PRELID1
Other Names	PX19; PRELI; SBBI12; CGI-106
Accession No.	Swiss-Prot#: Q9Y255 Gene ID: 27166
Uniprot	Q9Y255
GeneID	27166;
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:25-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human prostate cancer tissue using #42978 at dilution 1/40,

## Background

This gene encodes a member of the late embryogenesis abundant motif-containing protein family. The encoded protein is localized to mitochondria and may function as a cytoprotectant by regulating cell death and differentiation. Alternative splicing results in multiple transcript variants encoding different isoforms. Several related pseudogenes have been identified.

---

Note: This product is for in vitro research use only