

# Recombinant Human CAMK4

Catalog No: #GP11310



Package Size: #GP11310-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CAMK4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 190-313 amino acids of human calcium/calmodulin-dependent protein kinase IV
Target Name	calcium/calmodulin-dependent protein kinase IV
Other Names	IV; caMK; CaMK IV; CaMK-GR
Accession No.	Swissprot:Q16566Gene Accession:BC016695
Uniprot	Q16566
GenID	814;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The product of this gene belongs to the serine/threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. This enzyme is a multifunctional serine/threonine protein kinase with limited tissue distribution, that has been implicated in transcriptional regulation in lymphocytes, neurons and male germ cells.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only