

# Recombinant Human NCAPD2

Catalog No: #GP11312



Package Size: #GP11312-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human NCAPD2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 300 amino acids of human non-SMC condensin I complex, subunit D2
Target Name	non-SMC condensin I complex, subunit D2
Other Names	CNAP1; CAP-D2; hCAP-D2
Accession No.	Swissprot:Q15021Gene Accession:BC028182
Uniprot	Q15021
GenID	9918;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Condensin complex subunit 1 (also known as chromosome-associated protein D2 (CAP-D2) or non-SMC condensin I complex subunit D2 (NCAPD2) or XCAP-D2 homolog) is a protein that in humans is encoded by the NCAPD2 gene. CAP-D2 is a subunit of condensin I, a large protein complex involved in chromosome condensation.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only