

NUDCD3 Antibody

Catalog No: #42980

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NUDCD3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NUDCD3 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human NUDCD3
Target Name	NUDCD3
Other Names	NudCL
Accession No.	Swiss-Prot#: Q8IVD9Gene ID: 23386
Uniprot	Q8IVD9
GeneID	23386;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colorectal cancer tissue using #42980 at dilution 1/200.



Immunohistochemical analysis of paraffin-embedded Human Esophagus cancer tissue using #42980 at dilution 1/200.

Background

The product of this gene functions to maintain the stability of dynein intermediate chain. Depletion of this gene product results in aggregation and degradation of dynein intermediate chain, mislocalization of the dynein complex from kinetochores, spindle microtubules, and spindle poles, and loss of gamma-tubulin from spindle poles. The protein localizes to the Golgi apparatus during interphase, and levels of the protein increase after the G1/S transition.

Note: This product is for in vitro research use only