

Recombinant Human BSG

Catalog No: #GP11331



Package Size: #GP11331-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BSG
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 31-209 amino acids of human basigin
Target Name	basigin
Other Names	M6; OK; 5F7; TCSF; CD147; EMMPRIN
Accession No.	Swissprot:P35613Gene Accession:BC009040
Uniprot	P35613
GeneID	682;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a plasma membrane protein that is important in spermatogenesis, embryo implantation, neural network formation, and tumor progression. The encoded protein is also a member of the immunoglobulin superfamily. Multiple transcript variants encoding different isoforms have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only