

## DNAAF5 Antibody

Catalog No: #42981

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	DNAAF5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DNAAF5 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human DNAAF5
Target Name	DNAAF5
Other Names	CILD18; HEATR2
Accession No.	Swiss-Prot#: Q86Y56 Gene ID: 54919
Uniprot	Q86Y56
GeneID	54919;
Concentration	1.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:25-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #42981 at dilution 1/35,

## Background

The protein encoded by this gene is essential for the preassembly or stability of axonemal dynein arms, and is found only in organisms with motile cilia and flagella. Mutations in this gene are associated with primary ciliary dyskinesia-18, a disorder characterized by abnormalities of motile cilia. Alternatively spliced transcript variants have been found for this gene.

---

Note: This product is for in vitro research use only