

# Recombinant Human CD244

Catalog No: #GP11346



Package Size: #GP11346-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CD244
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 251-365 amino acids of human CD244 molecule, natural killer cell receptor 2B4
Target Name	CD244 molecule, natural killer cell receptor 2B4
Other Names	2B4; NAIL; Nmrk; NKR2B4; SLAMF4
Accession No.	Swissprot:Q9BZW8Gene Accession:BC028073
Uniprot	Q9BZW8
GenID	51744;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a cell surface receptor expressed on natural killer (NK) cells (and some T cells) that mediate non-major histocompatibility complex (MHC) restricted killing. The interaction between NK-cell and target cells via this receptor is thought to modulate NK-cell cytolytic activity. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only