

Recombinant Human BTLA

Catalog No: #GP11348



Package Size: #GP11348-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BTLA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 31-157 amino acids of human BTLA
Target Name	B and T lymphocyte associated
Other Names	BTLA1; CD272
Accession No.	Swissprot:Q7Z6A9Gene Accession:BC107091
Uniprot	Q7Z6A9
GeneID	151888;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the immunoglobulin superfamily. The encoded protein contains a single immunoglobulin (Ig) domain and is a receptor that relays inhibitory signals to suppress the immune response. Alternative splicing results in multiple transcript variants. Polymorphisms in this gene have been associated with an increased risk of rheumatoid arthritis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only