

Recombinant Human CD40LG

Catalog No: #GP11352



Package Size: #GP11352-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CD40LG
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 47-261 amino acids of human CD40 ligand
Target Name	CD40 ligand
Other Names	IGM; IMD3; TRAP; gp39; CD154; CD40L; HIGM1; T-BAM; TNFSF5; hCD40L
Accession No.	Swissprot:P29965Gene Accession:BC071754
Uniprot	P29965
GeneID	959;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is expressed on the surface of T cells. It regulates B cell function by engaging CD40 on the B cell surface. A defect in this gene results in an inability to undergo immunoglobulin class switch and is associated with hyper-IgM syndrome

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only