

# Recombinant Human CDC37L1

Catalog No: #GP11356



Package Size: #GP11356-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CDC37L1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	cell division cycle 37-like 1
Other Names	HARC; CDC37B; RP11-6J24.5
Accession No.	Swissprot:Q7L3B6Gene Accession:BC014133
Uniprot	Q7L3B6
GeneID	55664;
Storage	-20~-80°C, pH 7.6 PBS

## Background

CDC37L1 is a cytoplasmic phosphoprotein that exists in complex with HSP90 (HSPCA; MIM 140571) as well as several other proteins involved in HSP90-mediated protein folding

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only