

Recombinant Human CDK5RAP3

Catalog No: #GP11358



Package Size: #GP11358-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CDK5RAP3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 248-506 amino acids of human CDK5 regulatory subunit associated protein 3
Target Name	CDK5 regulatory subunit associated protein 3
Other Names	C53; IC53; LZAP; HSF-27; MST016; PP1553; OK/SW-cl.114
Accession No.	Swissprot:Q96JB5Gene Accession:BC009957
Uniprot	Q96JB5
GenID	80279;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that has been reported to function in signaling pathways governing transcriptional regulation and cell cycle progression. It may play a role in tumorigenesis and metastasis. A pseudogene of this gene is located on the long arm of chromosome 20. Alternative splicing results in multiple transcript variants that encode different isoforms.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only