

Recombinant Human CENPO

Catalog No: #GP11364



Package Size: #GP11364-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CENPO
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	centromere protein O
Other Names	CENP-O; MCM21R; ICEN-36
Accession No.	Swissprot:Q9BU64Gene Accession:BC002870
Uniprot	Q9BU64
GeneID	79172;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a component of the interphase centromere complex. The encoded protein is localized to the centromere throughout the cell cycle and is required for bipolar spindle assembly, chromosome segregation and checkpoint signaling during mitosis. Alternatively spliced transcript variants encoding multiple protein isoforms have been observed for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only