

## ZNF22 Antibody

Catalog No: #42986

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF22 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF22 protein.
Immunogen Description	Fusion protein of human ZNF22
Target Name	ZNF22
Other Names	KOX15; HKR-T1; ZNF422; Zfp422
Accession No.	Swiss-Prot#: P17026Gene ID: 7570
Uniprot	P17026
GeneID	7570;
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Lung cancer tissue using #42986 at dilution 1/200.



Immunohistochemical analysis of paraffin-embedded Human Thyroid cancer tissue using #42986 at dilution 1/200.

## Background

Binds DNA through the consensus sequence 5'-CAATG-3'. May be involved in transcriptional regulation and may play a role in tooth formation (By similarity).

Note: This product is for in vitro research use only