

# Recombinant Human CNTROB

Catalog No: #GP11389



Package Size: #GP11389-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CNTROB
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human centrobins, centrosomal BRCA2 interacting protein
Target Name	centrobins, centrosomal BRCA2 interacting protein
Other Names	LIP8; PP1221
Accession No.	Swissprot:Q8N137Gene Accession:BC021134
Uniprot	Q8N137
GeneID	116840;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a centrosomal protein that interacts with BRCA2, and is required for centriole duplication and cytokinesis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only