

## Recombinant Human CPA2

Catalog No: #GP11395



Package Size: #GP11395-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CPA2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 178-350 amino acids of human carboxypeptidase A2 (pancreatic)
Target Name	carboxypeptidase A2 (pancreatic)
Accession No.	Swissprot:P48052Gene Accession:BC007009
Uniprot	P48052
GeneID	1358;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Three different forms of human pancreatic procarboxypeptidase A have been isolated. The encoded protein represents the A2 form, which is a monomeric protein with different biochemical properties from the A1 and A3 forms. The A2 form of pancreatic procarboxypeptidase acts on aromatic C-terminal residues and is a secreted protein.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only