

# Recombinant Human CPB1

Catalog No: #GP11396



Package Size: #GP11396-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CPB1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 111-352 amino acids of human carboxypeptidase B1 (tissue)
Target Name	carboxypeptidase B1 (tissue)
Other Names	CPB; PASP; PCPB
Accession No.	Swissprot:P15086Gene Accession:BC015338
Uniprot	P15086
GenID	1360;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Three different procarboxypeptidases A and two different procarboxypeptidases B have been isolated. The B1 and B2 forms differ from each other mainly in isoelectric point. Carboxypeptidase B1 is a highly tissue-specific protein and is a useful serum marker for acute pancreatitis and dysfunction of pancreatic transplants. It is not elevated in pancreatic carcinoma.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only