

## FAM3D Antibody

Catalog No: #42987

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FAM3D Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FAM3D protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human FAM3D
Target Name	FAM3D
Other Names	EF7; OIT1
Accession No.	Swiss-Prot#: Q96BQ1 Gene ID: 131177
Uniprot	Q96BQ1
GeneID	131177;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Colorectal cancer tissue using #42987 at dilution 1/200,



Immunohistochemical analysis of paraffin-embedded Human Lung cancer tissue using #42987 at dilution 1/200,

## Background

FAM3D, also known as EF-7, belongs to the novel cytokine-like family whose members are secreted four-helix bundle proteins. Human FAM3D is abundantly expressed in placenta as a 224 amino acid (aa) residue precursor with a 25 aa signal sequence. The function of this protein remains unknown.

Note: This product is for in vitro research use only