

Recombinant Human CSNK1D

Catalog No: #GP11405



Package Size: #GP11405-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CSNK1D
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 250 amino acids of human septin 7
Target Name	casein kinase 1, delta
Other Names	ASPS; HCKID; FASPS2; CKIdelta
Accession No.	Swissprot:P48730Gene Accession:BC003558
Uniprot	P48730
GeneID	1453;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is a member of the casein kinase I (CKI) gene family whose members have been implicated in the control of cytoplasmic and nuclear processes, including DNA replication and repair. The encoded protein is highly similar to the mouse and rat CK1 delta homologs. Two transcript variants encoding different isoforms have been found for this gene

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only