

Recombinant Human CTAGE5

Catalog No: #GP11406



Package Size: #GP11406-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CTAGE5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 321-501 amino acids of human CTAGE family, member 5
Target Name	CTAGE family, member 5
Other Names	MEA6; MGEA; MGEA6; MGEA11
Accession No.	Swissprot:O15320Gene Accession:BC064355
Uniprot	O15320
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a tumor-associated antigen found in cutaneous T-cell lymphoma and several other cancers. Autoantibodies against the encoded protein have been found in some cancers. Several transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only