

## Recombinant Human RBBP8

Catalog No: #GP11407



Package Size: #GP11407-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RBBP8
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 603-902 amino acids of human RBBP8
Target Name	retinoblastoma binding protein 8
Other Names	RIM; COM1; CTIP; JWDS; SAE2; SCKL2
Accession No.	Swissprot:Q99708Gene Accession:BC030590
Uniprot	Q99708
GeneID	5932;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is a ubiquitously expressed nuclear protein. It is found among several proteins that bind directly to retinoblastoma protein, which regulates cell proliferation. This protein complexes with transcriptional co-repressor CTBP. It is also associated with BRCA1 and is thought to modulate the functions of BRCA1 in transcriptional regulation, DNA repair, and/or cell cycle checkpoint control. It is suggested that this gene may itself be a tumor suppressor acting in the same pathway as BRCA1.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only