

Recombinant Human CTNNBIP1

Catalog No: #GP11408

Package Size: #GP11408-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CTNNBIP1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	catenin, beta interacting protein 1
Other Names	ICAT
Accession No.	Swissprot:Q9NSA3Gene Accession:BC014300
Uniprot	Q9NSA3
GeneID	56998;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene binds CTNNB1 and prevents interaction between CTNNB1 and TCF family members. The encoded protein is a negative regulator of the Wnt signaling pathway. Two transcript variants encoding the same protein have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only