

## Recombinant Human CYP2W1

Catalog No: #GP11418



Package Size: #GP11418-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CYP2W1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from C terminal 250 amino acids of human cytochrome P450, family 2, subfamily W, polypeptide 1
Target Name	cytochrome P450, family 2, subfamily W, polypeptide 1
Accession No.	Swissprot:Q8TAV3Gene Accession:BC025761
Uniprot	Q8TAV3
GeneID	54905;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only