

Recombinant Human HSPA4

Catalog No: #GP11426



Package Size: #GP11426-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HSPA4
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	heat shock 70kDa protein 4
Other Names	RY; APG-2; HSPH2; hsp70; hsp70RY; HS24/P52
Accession No.	Swissprot:P34932Gene Accession:BC002526
Uniprot	P34932
GeneID	3308;
Storage	-20~-80°C, pH 7.6 PBS

Background

The heat shock proteins (HSPs) comprise a group of highly conserved, abundantly expressed proteins with diverse functions, which include the assembly and sequestering of multiprotein complexes, transportation of nascent polypeptide chains across cellular membranes and regulation of protein folding. Heat shock proteins (also known as molecular chaperones) fall into six general families: HSP 90, HSP 70, HSP 60, the low molecular weight HSPs, the immunophilins and the HSP 110 family. The HSP 110 family (also known as the HSP 105 family) is composed of HSP 105, Apg-1 and Apg-2.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only