

Recombinant Human CYR61

Catalog No: #GP11429



Package Size: #GP11429-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CYR61
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 25-381 amino acids of human cysteine-rich, angiogenic inducer, 61
Target Name	cysteine-rich, angiogenic inducer, 61
Other Names	CCN1; GIG1; IGFBP10
Accession No.	Swissprot:O00622Gene Accession:BC001271
Uniprot	O00622
GenID	3491;
Storage	-20~-80°C, pH 7.6 PBS

Background

The secreted protein encoded by this gene is growth factor-inducible and promotes the adhesion of endothelial cells. The encoded protein interacts with several integrins and with heparan sulfate proteoglycan. This protein also plays a role in cell proliferation, differentiation, angiogenesis, apoptosis, and extracellular matrix formation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only