

Recombinant Human BRINP1

Catalog No: #GP11439



Package Size: #GP11439-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BRINP1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	bone morphogenetic protein/retinoic acid inducible neural-specific 1
Other Names	DBC1; FAM5A; DBCCR1
Accession No.	Swissprot:O60477Gene Accession:BC021560
Uniprot	O60477
GeneID	1620;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene is located within a chromosomal region that shows loss of heterozygosity in some bladder cancers. It contains a 5' CpG island that may be a frequent target of hypermethylation, and it may undergo hypermethylation-based silencing in some bladder cancers.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only