

Recombinant Human DCTN5

Catalog No: #GP11448



Package Size: #GP11448-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DCTN5
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	dynactin 5 (p25)
Accession No.	Swissprot:Q9BTE1 Gene Accession:BC004191
Uniprot	Q9BTE1
GeneID	84516;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a subunit of dynactin, a component of the cytoplasmic dynein motor machinery involved in minus-end-directed transport. The encoded protein is a component of the pointed-end subcomplex and is thought to bind membranous cargo. A pseudogene of this gene is located on the long arm of chromosome 1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only