

## RNF141 Antibody

Catalog No: #42992

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	RNF141 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RNF141 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein of human RNF141
Target Name	RNF141
Other Names	ZFP26; ZNF230
Accession No.	Swiss-Prot#: Q8WVD5Gene ID: 50862
Uniprot	Q8WVD5
GeneID	50862;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:25-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human prostate cancer tissue using #42992 at dilution 1/25.



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #42992 at dilution 1/25.

## Background

The protein encoded by this gene contains a RING finger, a motif known to be involved in protein-DNA and protein-protein interactions. Abundant expression of this gene was found in the testicular tissue of fertile men, but was not detected in azoospermic patients. Studies of the mouse counterpart suggest that this gene may function as a testis specific transcription factor during spermatogenesis.

Note: This product is for in vitro research use only